

Chip-scale mid-infrared spectroscopy using electrically-pumped frequency comb sources

Lukasz A. Sterczewski^{1,2}, Tzu-Ling Chen³, Douglas C. Ober^{2,3}, Charles R. Markus³, Chadwick L. Canedy⁴, Igor Vurgaftman⁴, Clifford Frez², Mathieu Fradet², Siamak Forouhar², Mitchio Okumura³, Jerry R. Meyer⁴, and Mahmood Bagheri²

¹Faculty of Electronics, Photonics and Microsystems, Wroclaw University of Science and Technology, Wroclaw 50370, Poland

²Jet Propulsion Laboratory, California Institute of Technology, Pasadena, CA 91109, USA

³Division of Chemistry and Chemical Engineering, California Institute of Technology, Pasadena, CA, 91125, USA

⁴Optical Sciences Division, Naval Research Laboratory, Washington, DC 20375, USA

lukasz.sterczewski@pwr.edu.pl

Abstract: Interband cascade- and diode laser optical frequency combs with gigahertz repetition rates are discussed as sources for mode-resolved, broadband mid-infrared spectroscopy. Their millimeter-scale size and sub-watt power consumption are appealing for use in portable instruments.

The mid-wave infrared spectral region (MWIR, 3–5 μm) has come to play a critical role in sensing environmentally-relevant trace gases. In particular, chip-scale broadband coherent emitters with optical frequency comb (OFC) properties are highly-desirable for these applications. For spectroscopy, the discrete pattern of sharp, equidistant MWIR lines enables the multi-heterodyne (dual-comb) technique (Fig. 1a) to resolve individual comb lines on a nanosecond time scale with high signal-to-noise ratio (SNR). Especially attractive is the emerging interband cascade laser (ICL) OFC platform [1], which combines THz-wide spectral coverage and sub-MHz optical linewidth with the potential for realizing a complete room-temperature MWIR spectroscopic system that integrates a GHz-frequency interband cascade detector (see Fig. 1a) processed from the same wafer material as the ICL [2].

Figure 1b illustrates that the narrow optical linewidth of an ICL OFC also makes it well-suited for cavity-enhanced sensing techniques. The Vernier effect has been used in a centimeter-scale, high-finesse cavity to simultaneously enhance the optical path to 30 m and selectively filter the individual comb lines for broadband, millisecond-scale sensing of hydrocarbons at a wavelength near 3.6 μm [3]. This approach obviates the requirements for a fast photodetector and matched second comb source, which become even more difficult to fulfil in such challenging spectral regions as the long-wave infrared (LWIR) and terahertz (THz).

We will also discuss OFCs based on type-I GaSb Fabry-Pérot diode lasers, which have been shown to generate self-starting OFCs around 3.05 μm, with GHz repetition rate and frequency-modulated (FM) emission characteristics [4]. Although the timing stability is lower than for ICL OFCs, high spectral flatness and broad spectral tunability are combined with sub-MHz to MHz optical linewidth. This class of battery-operated OFC sources is well suited to the power-constrained detection of space-relevant molecules such as hydrogen cyanide (HCN) and acetylene (C₂H₂).




Fig. 1. Two topologies for chip-based, mode-resolved MWIR OFC spectroscopy: (a) Dual-comb spectrometer with a mercury-cadmium-telluride (MCT) or interband cascade (IC) detector. (b) Vernier spectrometer with cavity enhancement.

Acknowledgement

This work was supported under and was in part performed at the Jet Propulsion Laboratory (JPL), California Institute of Technology, under contract with NASA. L. A. Sterczewski's research was supported by an appointment to the NASA Postdoctoral Program at JPL, administered by Universities Space Research Association under contract with NASA. L. A. Sterczewski also acknowledges funding from the European Union's Horizon 2020 research and innovation programme under the Marie Skłodowska-Curie grant agreement No 101027721. The NRL authors acknowledge support from the Office of Naval Research.

References

- [1] M. Bagheri *et al.*, "Passively mode-locked interband cascade optical frequency combs," *Sci. Rep.* **8**, 3322 (2018).
- [2] L. A. Sterczewski *et al.*, "Mid-Infrared Dual-Comb Spectroscopy with Room-Temperature Bi-Functional Interband Cascade Lasers and Detectors," *Appl. Phys. Lett.* **116**, 141102 (2020).
- [3] L. A. Sterczewski *et al.*, "Cavity-Enhanced Vernier Spectroscopy with a Chip-Scale Mid-Infrared Frequency Comb," *ACS Photonics* (in press), (2022).
- [4] L. A. Sterczewski *et al.*, "Battery-operated mid-infrared diode laser frequency combs," *ArXiv211203964 Phys.* (2021).